



UniversityHospital Heidelberg

Ramin.Khoramnia@med.uni-heidelberg.de www.ivcrc.com www.djapplelab.com

Analysis of Opacified Hydrophilic IOLs After DSAEK

Khoramnia R.¹, De Cock R.², Emesz M.³, Auffarth G.U.¹

¹David J Apple Laboratory, Department of Ophthalmology, University of Heidelberg, Heidelberg, Germany (Chairman: G.U. Auffarth, MD, PhD, FEBO)

²Kent & Canterbury Hospital, Canterbury, UK

³Department of Ophthalmology, Paracelsus Medical University, Salzburg, Austria

IVCRC / DJ Apple Laboratory were supported by ...

Financial disclosure: 2013/2014

Alcon^{1,2,3}

Allergan^{1,4}

AMO^{1,2,3,4}

Bausch & Lomb / Technolas Perfect Vision^{1,2,3,4}

Bayer²

Carl Zeiss Meditec^{1,2}

Dr. Schmidt Intraocularlinsen^{1,3}

Geuder²

Novartis^{1,2}

Oculentis^{1,2,3}

Ophtec²

Physiol^{1,2}

Powervision¹

Rayner^{1,2,3}



1 = Research grant; 2 = Travel grant; 3 = Lecture honorarium; 4 = Consulting

Purpose

- Calcification of hydrophilic IOLs is a rare complication.¹
- The injection of gas into the anterior chamber during DSAEK seems to increase the risk for opacifications.²⁻⁴
- Granular deposits below the anterior surface of the IOLs can be responsible for such a decrease in visual acuity that patients require IOL exchange.
- The composition of the deposits can be analyzed by means of x-ray spectroscopy.¹

- 1.) Khoramnia R, Salgado JP, Auffarth GU, Schmidt S, Wegner A, Kobuch KA, Winkler von Mohrenfels C. Opacification of a hydrophilic intraocular lens 4 years after cataract surgery. A biomaterial analysis. *Ophthalmologe*. 2012 May;109(5):483-6.
- 2.) Dhital A, Spalton DJ, Goyal S, Werner L. Calcification in hydrophilic intraocular lenses associated with injection of intraocular gas. *Am J Ophthalmol*. 2012 Jun;153(6):1154-60.e1.
- 3.) Werner L, Wilbanks G, Ollerton A, Michelson J. Localized calcification of hydrophilic acrylic intraocular lenses in association with intracameral injection of gas. *J Cataract Refract Surg*. 2012 Apr;38(4):720-1.
- 4.) Schmidinger G, Pemp B, Werner L. Opacification of an intraocular lens : Calcification of hydrophilic intraocular lenses after gas tamponade of the anterior chamber. *Ophthalmologe*. 2013 Nov;110(11):1066-8.




Methods

- Six hydrophilic IOLs with opacification after DSAEK have been evaluated in this ongoing analysis so far.
- The explanted IOLs were cut in half.
- One half was stained (Alizarin red and von Kossa) and examined using light microscopy.
- The other half was examined using scanning electron microscopy.
- The composition of the deposits was analyzed by means of x-ray spectroscopy.

Results: Light microscopy

IOL 1



Overview, Alizarin red



Alizarin red, orig. magnif. x 4



Alizarin red, orig. magnif. x 10



von Kossa, orig. magnif. x 40

IOL 2



Overview



Alizarin red, orig. magnif. x 4




Alizarin red, orig. magnif. x 40



von Kossa, orig. magnif. x 40

IOL 3



Overview



Unstained, orig. magnif. x 40



Alizarin red, orig. magnif. x 4



von Kossa, orig. magnif. x 40

IOL 4



Overview



Alizarin red, orig. magnif. x 10




Alizarin red, orig. magnif. x 40



von Kossa, orig. magnif. x 40


IOL 5



Overview



Alizarin red, orig. magnif. x 4



Alizarin red, orig. magnif. x 20



von Kossa, orig. magnif. x 20

IOL 6



Overview



Unstained, orig. magnif. x 4




Alizarin red, orig. magnif. x 4




von Kossa, orig. magnif. x 20


Results: Scanning Electron Microscopy




IOL 1




IOL 2




IOL 3



IOL 4




IOL 5




IOL 6

Numerous fine, granular deposits within the optic of the IOL, which were distributed in a line parallel to the anterior surface of the IOL.

Results: Element Mapping




IOL 4



IOL 5

Results: X-Ray Spectroscopy

IOL 5



Calciumphosphate

Element	Series	unn. C [wt. %]	norm. C [wt. %]	Atom. C [at. %]	Error (1 Sigma) [wt. %]
Carbon	K-series	58.55	58.55	68.25	7.00
Oxygen	K-series	32.28	32.28	28.25	4.34
Phosphorus	K-series	2.38	2.38	1.08	0.12
Calcium	K-series	6.53	6.53	2.28	0.22
Chlorine	K-series	0.04	0.04	0.02	0.03
Sodium	K-series	0.22	0.22	0.14	0.05
Total:		100.00	100.00	100.00	

Discussion

- The cause for calcification of IOLs has not yet been determined.
- Opacifications are often localized in the pupillary area of the IOLs.
- The surface of the IOLs might be affected by the gas bubble
→ development of a crystallization nucleus¹
- The altered blood-aqueous barrier might co-act to trigger secondary IOL calcification.²

1.) Schmidinger G, Pemp B, Werner L. Opacification of an intraocular lens : Calcification of hydrophilic intraocular lenses after gas tamponade of the anterior chamber. *Ophthalmologe*. 2013 Nov;110(11):1066-8.

2.) Dhital A, Spalton DJ, Goyal S, Werner L. Calcification in hydrophilic intraocular lenses associated with injection of intraocular gas. *Am J Ophthalmol*. 2012 Jun;153(6):1154-60.e1.

Summary

- The injection of gas into the anterior chamber during DSAEK seems to increase the risk for opacifications.
- Calciumphosphate can be stained using Alizarin red or von Kossa.
- The exact chemical composition of the deposits can be assessed using X-ray spectroscopy.

Summary

- The incidence of opacifications of hydrophilic IOLs after DSAEK does not seem to be very high.
- Nevertheless, the implantation of a hydrophilic IOLs does not seem to be the ideal solution in patients who might require corneal surgery later on.¹
- An IOL exchange and the associated risks could be avoided in such patients by use of hydrophobic IOLs.¹

1.) Schmidinger G, Pemp B, Werner L. Opacification of an intraocular lens : Calcification of hydrophilic intraocular lenses after gas tamponade of the anterior chamber. *Ophthalmologe*. 2013 Nov;110(11):1066-8.



Zertifiziert
DIN EN ISO
9001:2008



G.U. Auffarth, MD, PhD, FEBO

M.P. Holzer, MD, FEBO
F.T.A. Kretz, MD, FEBO
K. Linz, MD

T.M. Rabsilber, MD
F.N. Auerbach, MD
A. Fitting, MSc

R. Khoramnia, MD, FEBO
M. Safwat Attia, MD
R. Willrich Amroussi, MA